

Product Datasheet

Anti-Cyclic Citrullinated Peptide (anti-CCP) Concentrate 194-75

| | |
|----------------------------|--|
| Product Name | Anti-Cyclic Citrullinated Peptide (anti-CCP) Concentrate |
| Catalog Number | 194-75 |
| Alternative Names | Anti-citrullinated protein antibody, aCCP, anti-CCP, CCP, Cyclic Citrullinated Peptide |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid |
| Concentration | Reported in mg/mL (Biuret) |
| Storage | -20°C |
| Purity | Concentrate |
| Source | Human Plasma |
| Biological Activity | $\geq > 5,000$ IU/mL (anti-CCP IgG ELISA) |