

## Product Datasheet

### Colipase 197-25

<b>Product Name</b>	Colipase
<b>Catalog Number</b>	197-25
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Source</b>	Porcine Pancreas
<b>Molecular Weight</b>	~10,300
<b>Biological Activity</b>	4,000-12,000 U/mg
<b>Specific Activity</b>	Reported in U/mg protein
<b>Unit Definition</b>	One unit of colipase corresponds to the amount of cofactor that increases bile salt-inhibited lipase activity by 1 enzyme unit