

Product Datasheet

Low Density Lipoprotein (LDL) - MDA Modified 360-27

Product Name Low Density Lipoprotein (LDL) - MDA Modified

Catalog Number 360-27

Alternative Names MDA modified LDL

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration ≥ 1 mg/mL (Coomassie)

Storage 2-8 °C

Purity $\geq 95\%$ (SDS-PAGE)

Source Human Plasma

Molecular Weight 2.75 million (ref. Tetz second Edition)