

Product Datasheet

Low Density Lipoprotein (LDL) 360-30

Product Name	Low Density Lipoprotein (LDL)
Catalog Number	360-30
Alternative Names	LDL, LDL-C, Beta-Lipoprotein, Concentrate, Fraction
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	5 - 8 mg/mL (Modified Lowry)
Storage	2-8 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Molecular Weight	2.75 million (ref. Tetz second Edition)
Triglyceride Ratio	≤ 0.50 (Triglycerides/Total Cholesterol)