

Product Datasheet

Low Density Lipoprotein (OxLDL) - Oxidized 360-31

Product Name Low Density Lipoprotein (OxLDL) - Oxidized

Catalog Number 360-31

Alternative Names OxLDL, Ox-LDL

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration 5 - 8 mg/mL (Modified Lowry)

Storage 2-8 °C

Purity ≥ 95% (SDS-PAGE before oxidation)

Source Human Plasma

Triglyceride Ratio ≤ 0.50 (Triglycerides/Total Cholesterol)