

## Product Datasheet

### Low Density Lipoprotein (OxLDL) - Oxidized 360-31

**Product Name** Low Density Lipoprotein (OxLDL) - Oxidized

**Catalog Number** 360-31

**Alternative Names** OxLDL, Ox-LDL

**Brand** Lee Biosolutions

**Form/Appearance** Liquid

**Concentration** 5 - 8 mg/mL (Modified Lowry)

**Storage** 2-8 °C

**Purity** ≥ 95% (SDS-PAGE before oxidation)

**Source** Human Plasma

**Triglyceride Ratio** ≤ 0.50 (Triglycerides/Total Cholesterol)