

Product Datasheet

High Density Lipoprotein (HDL) 361-25

Product Name	High Density Lipoprotein (HDL)
Catalog Number	361-25
Alternative Names	HDL, HDL-C, HDL Concentrate, HDL Fraction, alpha-Lipoprotein
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1.0 - 5.0 mg/mL (Lowry)
Storage	2-8 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Molecular Weight	175,000-360,000 (Tetz second Edition)