

Product Datasheet

Rheumatoid Factor (RF) Concentrate ? 5 kIU/mL 508-10

Product Name	Rheumatoid Factor (RF) Concentrate \geq 5 kIU/mL
Catalog Number	508-10
Alternative Names	RF, RF Concentrate
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	Reported in mg/mL (Biuret)
Storage	-20°C
Purity	Concentrate
Source	Human Defibrinated Plasma (multiple, pooled donors)
Biological Activity	\geq 5,000 IU/mL (Turbidimetric Immunoassay)
Specific Activity	Reported in U/mg protein (Biuret)