

## Product Datasheet

### Rheumatoid Factor (RF) Concentrate ? 5 kIU/mL 508-10

<b>Product Name</b>	Rheumatoid Factor (RF) Concentrate $\geq$ 5 kIU/mL
<b>Catalog Number</b>	508-10
<b>Alternative Names</b>	RF, RF Concentrate
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Reported in mg/mL (Biuret)
<b>Storage</b>	-20°C
<b>Purity</b>	Concentrate
<b>Source</b>	Human Defibrinated Plasma (multiple, pooled donors)
<b>Biological Activity</b>	$\geq$ 5,000 IU/mL (Turbidimetric Immunoassay)
<b>Specific Activity</b>	Reported in U/mg protein (Biuret)