

## Product Datasheet

### Recombinant IGF-1 antigen 610173

<b>Product Name</b>	Recombinant IGF-1 antigen
<b>Catalog Number</b>	610173
<b>Description</b>	Recombinant human IGF-1 antigen with a C-terminal histidine tag. Chromatographically purified.
<b>Species/Host</b>	E. coli
<b>Alternative Names</b>	Insulin-like Growth Factor 1, Human Insulin-like Growth Factor 1
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Buffer</b>	50 mM trisodium citrate, 0.9 % NaCl, pH 6.0
<b>Molecular Weight</b>	Predicted molecular weight: 9 kDa
<b>Shelf Life</b>	Unspecified