

Product Datasheet

Recombinant IGF-1 antigen 610173

Product Name	Recombinant IGF-1 antigen
Catalog Number	610173
Description	Recombinant human IGF-1 antigen with a C-terminal histidine tag. Chromatographically purified.
Species/Host	E. coli
Alternative Names	Insulin-like Growth Factor 1, Human Insulin-like Growth Factor 1
Brand	Medix Biochemica
Form/Appearance	Liquid
Storage	-20°C
Buffer	50 mM trisodium citrate, 0.9 % NaCl, pH 6.0
Molecular Weight	Predicted molecular weight: 9 kDa
Shelf Life	Unspecified