

## Product Datasheet

### Recombinant IGF-1 antigen 610173

|                          |  |
|--------------------------|--|
| <b>Product Name</b>      | Recombinant IGF-1 antigen  |
| <b>Catalog Number</b>    | 610173   |
| <b>Description</b>       | Recombinant human IGF-1 antigen with a C-terminal histidine tag. Chromatographically purified. |
| <b>Species/Host</b>      | E. coli  |
| <b>Alternative Names</b> | Insulin-like Growth Factor 1, Human Insulin-like Growth Factor 1                               |
| <b>Brand</b>             | Medix Biochemica   |
| <b>Form/Appearance</b>   | Liquid   |
| <b>Storage</b>           | -20°C  |
| <b>Buffer</b>            | 50 mM trisodium citrate, 0.9 % NaCl, pH 6.0  |
| <b>Molecular Weight</b>  | Predicted molecular weight: 9 kDa  |
| <b>Shelf Life</b>        | Unspecified  |