

## Product Datasheet

### Carbapenemase (KPC) Antibody HM1441

<b>Product Name</b>	Carbapenemase (KPC) Antibody
<b>Catalog Number</b>	HM1441
<b>Species/Host</b>	Mouse
<b>Alternative Names</b>	KPC
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	HM1440 or HM1442
<b>Preservatives</b>	0.09% Sodium Azide