

Product Datasheet

Carbapenemase (KPC) Antibody HM1442

| | |
|--------------------------|---|
| Product Name | Carbapenemase (KPC) Antibody |
| Catalog Number | HM1442 |
| Species/Host | Mouse |
| Alternative Names | KPC |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | Phosphate Buffered Saline, pH 7.4 |
| Matched Pairs | HM1441 |
| Preservatives | 0.09% Sodium Azide |