

Product Datasheet

Carbapenemase (KPC) Antibody HM1442

Product Name	Carbapenemase (KPC) Antibody
Catalog Number	HM1442
Species/Host	Mouse
Alternative Names	KPC
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM1441
Preservatives	0.09% Sodium Azide