

Product Datasheet

Carbapenemase (NDM) Antibody HM1451

Product Name	Carbapenemase (NDM) Antibody
Catalog Number	HM1451
Species/Host	Mouse
Alternative Names	NDM
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM1452
Preservatives	0.09% Sodium Azide