

Product Datasheet

Treponema pallidum (TpN17) Antigen, Recombinant >90% LA813

Product Name	Treponema pallidum (TpN17) Antigen, Recombinant >90%
Catalog Number	LA813
Tested Applications	ELISA, CMIA
Alternative Names	Syphilis
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>90%
Source	
Buffer	25 mM Tris-HCl, pH 7.0 with 30 mM NaCl & 3 mM EDTA
Specificity	TpN17
Expression System	E. coli