

Product Datasheet

Cortisol-3-BSA Antigen LA790

Product Name	Cortisol-3-BSA Antigen
Catalog Number	LA790
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM1255 or HM1256
Preservatives	0.05% Sodium Azide