

Product Datasheet

Anti-KLH Azide Free Clone B-K60 857.460.000, 857.460.005

Product Name	Anti-KLH Azide Free Clone B-K60
Catalog Number	857.460.000, 857.460.005
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Protein chimeric KLH
Alternative Names	Keyhole limpet hemocyanin
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	IgG1
Clonality	Monoclonal
Specificity	KLH (Keyhole limpet hemocyanin)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells