

Product Datasheet

Human IL-2 EasySplit ELISpot Kit 856.001.001PCS

Product Name	Human IL-2 EasySplit ELISpot Kit
Catalog Number	856.001.001PCS, 856.001.002PCS, 856.001.005PCS
Description	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Brand	Diaclone
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-2.
Specificity	Recognizes natural human IL-2
References	Junwei, W. et al., Cancer Cell Int.,2016; 17: 10. https://www.ncbi.nlm.nih.gov/pubmed/27330408