

## Product Datasheet

### Human IL-2 EasySplit ELISpot Kit 856.001.001PCS

<b>Product Name</b>	Human IL-2 EasySplit ELISpot Kit
<b>Catalog Number</b>	856.001.001PCS, 856.001.002PCS, 856.001.005PCS
<b>Description</b>	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
<b>Tested Applications</b>	ELISpot
<b>Brand</b>	Diaclone
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with simian IL-2.
<b>Specificity</b>	Recognizes natural human IL-2
<b>References</b>	Junwei, W. et al., Cancer Cell Int.,2016; 17: 10. <a href="https://www.ncbi.nlm.nih.gov/pubmed/27330408">https://www.ncbi.nlm.nih.gov/pubmed/27330408</a>