

Product Datasheet

Human Perforin EasySplit ELISpot Kit 856.141.001PCS

Product Name	Human Perforin EasySplit ELISpot Kit
Catalog Number	856.141.001PCS, 856.141.002PCS, 856.141.005PCS
Description	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Brand	Diaclone
Specificity	Recognizes natural human Perforin
References	Van Gulck, E. R. et al, J. Virol,2008; 82(7): 3561-3573. https://www.ncbi.nlm.nih.gov/pubmed/18234800