

## Product Datasheet

### Human Perforin EasySplit ELISpot Kit 856.141.001PCS

<b>Product Name</b>	Human Perforin EasySplit ELISpot Kit
<b>Catalog Number</b>	856.141.001PCS, 856.141.002PCS, 856.141.005PCS
<b>Description</b>	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
<b>Tested Applications</b>	ELISpot
<b>Brand</b>	Diaclone
<b>Specificity</b>	Recognizes natural human Perforin
<b>References</b>	Van Gulck, E. R. et al, J. Virol,2008; 82(7): 3561-3573. <a href="https://www.ncbi.nlm.nih.gov/pubmed/18234800">https://www.ncbi.nlm.nih.gov/pubmed/18234800</a>