

Product Datasheet

Human TNF- α EasySplit ELISpot Kit 856.041.001PCS, 856.041.002PCS, 856.041.005PCS

Product Name	Human TNF- α EasySplit ELISpot Kit
Catalog Number	856.041.001PCS, 856.041.002PCS, 856.041.005PCS
Description	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Brand	Diaclone
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with macac TNF- α . No interference with TNFR I and TNFR II
Specificity	Recognizes natural human TNF- α
References	Godard B. et al., Hum. Immunol., 2004; 65(11): 1307-18 https://www.ncbi.nlm.nih.gov/pubmed/15556681 Hudak S. et al., J. Immuno., 2002; 169(3): 1189 - 1196 https://www.ncbi.nlm.nih.gov/pubmed/12133939