

Product Datasheet

Human IL-6 EasySplit ELISpot Kit 856.021.001PCS, 856.021.002PCS, 856.021.005PCS

Product Name	Human IL-6 EasySplit ELISpot Kit
Catalog Number	856.021.001PCS, 856.021.002PCS, 856.021.005PCS
Description	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Brand	Diaclone
Cross Reactivity	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Specificity	Recognizes natural human IL-6