

Product Datasheet

Murine IL-2 EasySplit ELISpot Kit **862.001.001PCS, 862.001.002PCS, 862.001.005PCS**

Product Name	Murine IL-2 EasySplit ELISpot Kit
Catalog Number	862.001.001PCS, 862.001.002PCS, 862.001.005PCS
Description	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Brand	Diaclone
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IL-2
References	Pollard, C. et al., Mol Ther.,2013; 21(1): 251-9 https://www.ncbi.nlm.nih.gov/pubmed/23011030