

---

## Product Datasheet

### **Rat IFN- $\gamma$ EasySplit ELISpot Kit** **867.011.001PCS, 867.011.002PCS, 867.011.005PCS**

<b>Product Name</b>	Rat IFN- $\gamma$ EasySplit ELISpot Kit
<b>Catalog Number</b>	867.011.001PCS, 867.011.002PCS, 867.011.005PCS
<b>Description</b>	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
<b>Tested Applications</b>	ELISpot
<b>Brand</b>	Diaclone
<b>Cross Reactivity</b>	No cross reactivity with other rat cytokines
<b>Specificity</b>	Recognizes natural rat IFN-g