

Product Datasheet

Rat IFN- γ EasySplit ELISpot Kit **867.011.001PCS, 867.011.002PCS, 867.011.005PCS**

Product Name	Rat IFN- γ EasySplit ELISpot Kit
Catalog Number	867.011.001PCS, 867.011.002PCS, 867.011.005PCS
Description	Diaclone EasySplit ELISpot kits include: Pre-coated strip plate Detection biotinylated antibody Streptavidin Alkaline Phosphatase Conjugate Bovine Serum Albumin Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Brand	Diaclone
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat IFN-g