

Product Datasheet

Anti-Fentanyl 13102 SPTN-5 100974

| | |
|----------------------------|--|
| Product Name | Anti-Fentanyl 13102 SPTN-5 |
| Catalog Number | 100974 |
| Tested Applications | FIA |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Storage | +2-8 °C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 95 % |
| Affinity constant | N/A (not applicable) |
| Associated Products | 170040 |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 6.4-6.9 |
| Cross Reactivity | Recognizes acetyl fentanyl Does not recognize norfentanyl Risperidone ≤ 1 % 9-hydroxyrisperidone ≤ 1 % |
| Specificity | Antibody recognizes fentanyl |
| Matched Pairs | Fentanyl-BSA Antigen 170040 (capture) |
| Shelf Life | Unspecified |