

Product Datasheet

Xylazine Antibody HM1436

Product Name	Xylazine Antibody
Catalog Number	HM1436
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Xylazine-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Xylazine, 3-hydroxy Xylazine and 4-hydroxy Xylazine 100%; 4-hydroxy xylazine-O-glucuronide 120%
Specificity	Xylazine and metabolites
Matched Pairs	LA815
Preservatives	0.05% Sodium Azide