

Product Datasheet

Xylazine Antibody HM1436

Product Name Xylazine Antibody

Catalog Number HM1436

Species/Host Mouse

Tested Applications LF, ELISA

Immunogen Xylazine-KLH

Brand EastCoast Bio

Form/Appearance Ion Exchange Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Buffer 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride

Cross Reactivity Xylazine, 3-hydroxy Xylazine and 4-hydroxy Xylazine 100%;
4-hydroxy xylazine-O-glucuronide 120%

Specificity Xylazine and metabolites

Matched Pairs LA815

Preservatives 0.05% Sodium Azide