

## Product Datasheet

### Xylazine-BSA Antigen LA815

**Product Name** Xylazine-BSA Antigen

**Catalog Number** LA815

**Tested Applications** LF, ELISA

**Form/Appearance** Liquid; Antigen-BSA conjugate

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Purity** ≥90% by SDS-PAGE

**Buffer** 10 mM Phosphate Buffer, pH 7.2

**Matched Pairs** HM1435, HM1436, HM1437

**Preservatives** 0.01% Sodium Azide