

## Product Datasheet

### Anti-EtG R12403 SPTNZ-1 140052

<b>Product Name</b>	Anti-EtG R12403 SPTNZ-1
<b>Catalog Number</b>	140052
<b>Tested Applications</b>	FIA
<b>Alternative Names</b>	Ethyl glucuronide, Ethanol
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	1.0 mg/ml (+/-10 %) (A280 nm, 1 mg/ml, 1 cm = 1.6)
<b>Storage</b>	2-8 °C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, with 0.05 % sulfobetaine and 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	7.0-7.8
<b>Cross Reactivity</b>	Not determined (N/D)
<b>Specificity</b>	Antibody recognizes ethyl glucuronide (EtG)
<b>Matched Pairs</b>	EtG-BSA antigen LA189
<b>Shelf Life</b>	Unspecified