

## Product Datasheet

### Anti-6-MAM R12703 SPTN-5 140041

<b>Product Name</b>	Anti-6-MAM R12703 SPTN-5
<b>Catalog Number</b>	140041
<b>Tested Applications</b>	FIA, Lateral Flow
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	Not Determined
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.4-6.6
<b>Cross Reactivity</b>	Morphine ≤ 0.01 %
<b>Specificity</b>	Antibody recognizes 6-monoacetylmorphine (6-MAM)
<b>Matched Pairs</b>	6-MAM-BSA Antigen 170031
<b>Shelf Life</b>	Unspecified