

Product Datasheet

Mycobacterium tuberculosis MPT64 (TB-MPT64) Antibody HM1423

Product Name	Mycobacterium tuberculosis MPT64 (TB-MPT64) Antibody
Catalog Number	HM1423
Tested Applications	LF
Immunogen	M. tuberculosis MPT64 (recombinant)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	≥95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM1424 (detection)
Preservatives	0.02% Sodium Azide