

## Product Datasheet

### Brain Natriuretic Peptide (BNP) Antigen ≥90% LA788

<b>Product Name</b>	Brain Natriuretic Peptide (BNP) Antigen ≥90%
<b>Catalog Number</b>	LA788
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized Antigen
<b>Storage</b>	Short Term (≤ 2 weeks): 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90%
<b>Source</b>	Synthetic
<b>Matched Pairs</b>	HM818 and HM819