

## Product Datasheet

### Human T-lymphotropic Virus 1 (HTLV-1 gp21) Antigen, Recombinant LA816

<b>Product Name</b>	Human T-lymphotropic Virus 1 (HTLV-1 gp21) Antigen, Recombinant
<b>Catalog Number</b>	LA816
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	IT
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	Short Term ( $\leq 1$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Epitope</b>	HTLV-1 gp21 (1-152 aa)
<b>Purity</b>	$\geq 90$ by SDS-PAGE
<b>Buffer</b>	Phosphate Buffered Saline
<b>Specificity</b>	
<b>Molecular Weight</b>	17.9 kDa