

Product Datasheet

Human T-lymphotropic Virus 1 (HTLV-1 gp46) Antigen, Recombinant LA817

Product Name	Human T-lymphotropic Virus 1 (HTLV-1 gp46) Antigen, Recombinant
Catalog Number	LA817
Species/Host	E. coli
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (≤ 1 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Epitope	HTLV-1 gp46 (21-313 aa)
Purity	≥ 90 by SDS-PAGE
Buffer	20 mM Phosphate Buffer, pH 7.5
Specificity	
Molecular Weight	34.2 kDa