

Product Datasheet

Aripiprazole Antibody HM1263

Product Name	Aripiprazole Antibody
Catalog Number	HM1263
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Aripiprazole-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C (short-term) -20°C (long-term)
Isotype	IgG1
Clonality	Monoclonal
Purity	≥90%
Buffer	0.015M Potassium Phosphate Buffer with 0.85% NaCl, pH 7.2
Cross Reactivity	Aripiprazole 100%; Clozapine, Olanzapine, Quetiapine, Risperidone, Ziprasidone, Amitriptyline, Clomipramine, Cyclobenzaprine, Desipramine, Doxepin, Imipramine, Maprotiline, Nordoxepin, Nortriptyline, Trimipramine <1%
Specificity	Aripiprazole
Matched Pairs	LA799 (capture)
Preservatives	0.05% Sodium Azide