

Product Datasheet

Methamphetamine Antibody HM1268

| | |
|----------------------------|---|
| Product Name | Methamphetamine Antibody |
| Catalog Number | HM1268 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | DEAE Purified Monoclonal Antibody |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | Amphetamine 0.8%, MDMA 0.3 % |
| Specificity | Methamphetamine |
| Matched Pairs | LA469 (capture) |
| Preservatives | 0.05% Sodium Azide |