

Product Datasheet

Methamphetamine Antibody HM1268

Product Name Methamphetamine Antibody

Catalog Number HM1268

Species/Host Mouse

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance DEAE Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Buffer 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride

Cross Reactivity Amphetamine 0.8%, MDMA 0.3 %

Specificity Methamphetamine

Matched Pairs LA469 (capture)

Preservatives 0.05% Sodium Azide