

Product Datasheet

Clozapine-BSA Antigen LA800

Product Name	Clozapine-BSA Antigen
Catalog Number	LA800
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM1272 or HM1273 (detection)
Preservatives	0.01% Sodium Azide