

Product Datasheet

Quetiapine Antibody HM1281

Product Name	Quetiapine Antibody
Catalog Number	HM1281
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Quetiapine-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥90%
Buffer	0.015M Potassium Phosphate Buffer with 0.85% NaCl, pH 7.2
Cross Reactivity	Quetiapine 100%; Aripiprazole, Clozapine, Olanzapine, Risperidone, Ziprasidone, Amitriptyline, Clomipramine, Cyclobenzaprine, Desipramine, Doxepin, Imipramine, Maprotiline, Nordoxepin, Nortriptyline, Trimipramine <1%
Specificity	Quetiapine
Matched Pairs	LA803 (capture)
Preservatives	0.05% Sodium Azide