

# Product Datasheet

## Serum, Reverse Triiodothyronine (T3) Samples 991-24-PS-RT3

Product Name	Serum, Reverse Triiodothyronine (T3) Samples
Catalog Number	991-24-PS-RT3
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.