

Product Datasheet

Serum, Triiodothyronine (T3) Uptake Samples 991-24-PS-UT3

Product Name Serum, Triiodothyronine (T3) Uptake Samples

Catalog Number 991-24-PS-UT3

Brand Lee Biosolutions

Form/Appearance Liquid

Storage -20°C

Source Single Human Donor

Shelf Life Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.