

## Product Datasheet

### Serum, Topiramate Samples 991-24-PS-TOPI

**Product Name** Serum, Topiramate Samples

**Catalog Number** 991-24-PS-TOPI

**Brand** Lee Biosolutions

**Form/Appearance** Liquid

**Storage** -20°C

**Source** Single Human Donor

**Shelf Life** Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.