

Product Datasheet

Coronavirus (COVID-19 NP; 2-419) Antigen, Recombinant >95% LA600

Product Name	Coronavirus (COVID-19 NP; 2-419) Antigen, Recombinant >95%
Catalog Number	LA600
Species/Host	E. coli
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>95%
Buffer	Phosphate Buffered Saline with 25 mM Potassium Carbonate, pH 9.0
Specificity	Wuhan-Hu-1
Molecular Weight	48 kDa
Matched Pairs	HM1054-HM1058, HM1063-HM1065, HM1143-HM1147
Expression System	E. coli
Preservatives	0.05% Sodium Azide