
Product Datasheet

Serum, Troponin I Samples 991-24-PS-TNI

Product Name Serum, Troponin I Samples

Catalog Number 991-24-PS-TNI

Species/Host Human

Brand Lee Biosolutions

Form/Appearance Liquid

Storage -20 °C

Source Single Human Donor

Shelf Life Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.