

Product Datasheet

Anti-Human CD141 Monoclonal Antibody, Azide Free Clone B-S24 853.850

Product Name Anti-Human CD141 Monoclonal Antibody, Azide Free Clone B-S24

Catalog Number 853.850

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA

Immunogen Eahy 926 cell lines

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Source Recombinant antibody from the original clone

Specificity Recognises thrombomodulin (TM), a 105 kDa protein