

Product Datasheet

Anti-RSV NP R13401 SPTN-5 140042

Product Name	Anti-RSV NP R13401 SPTN-5
Catalog Number	140042
Description	Recombinant rabbit antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.
Tested Applications	LF, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	Rabbit IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Native RSV antigen FC517
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-6.9
Cross Reactivity	Does not recognize Influenza A NP, Influenza B NP and SARS-CoV-2 NP.
Specificity	Antibody recognizes respiratory syncytial virus NP, strains A & B
Matched Pairs	FIA: R13403 (capture) - R13404 (detection) LF: R13404 (membrane) - V0671 (particles); R13401 - V0671
Shelf Life	Unspecified